



Topical Silicone Gels

Your physician has cleared you to begin using Silagen® silicone gel on your post-op scars. Clinical studies of silicone gel demonstrate that they are able to **flatten and soften scars, as well as improve redness, itching and pain**. They have also been shown to help **speed up the improvement in appearance** compared to untreated scars. Following the directions below **exactly as directed** will help you optimize your scar treatment for the best results.

You may begin using the gel once the wound has closed—no bleeding or oozing—and after any removable stitches have been removed or as directed by your physician

DIRECTIONS FOR USE

1. Ensure scar area is clean and dry, as well as free from creams, lotions, ointments and powders.
2. Apply a very small amount of gel to scar and gently massage.
3. Allow to dry.
4. Once dry, cosmetics and compression garments may be used over the scar.
5. Use twice daily for 12 weeks for optimal results. Patients who use Silagen consistently twice a day for 3 months get the best results.
6. If you are prone to forming excessive scars (hypertrophic or keloid) you may want to continue treatment for up to 6 months.
7. Do not use the gel under Silagen silicone sheets, strips or shapes. They will not stick to the skin if you do.

CONTACT OUR OFFICE WITH ANY QUESTIONS:

CLINICAL EVIDENCE
OF
IMPROVED
SCAR
APPEARANCE
WITH USE OF
SILAGEN
SILICONE GEL



BEFORE:
3 days after
Mohs surgery



AFTER:
Scar after 2 months of
using Silagen Silicone
Gel twice a day.



AFTER:
After 4 months of
using Silagen Silicone
Gel twice a day.